

	Screen		Confirm	
Drug	Methodology	Cutoff	GC-MS	LC-MS/MS
<b>Alcohol (Ethanol)</b>	EIA	0.025 gm/dL	.025 gm/dL (GC-FID)	
<b>Amphetamines</b>				
- Amphetamine				15 ng/mL
- Methamphetamine				15 ng/mL
- MDA				15 ng/mL
- MDMA				15 ng/mL
- MDEA				15 ng/mL
<b>Barbiturates</b>	EIA	50 ng/mL		
- Butalbital				25 ng/mL
- Pentobarbital				25 ng/mL
- Phenobarbital				25 ng/mL
- Secobarbital				25 ng/mL
<b>Benzodiazepines</b>	EIA	50 ng/mL		
- Alprazolam				0.5 ng/mL
- Chlordiazepoxide				0.5 ng/mL
- Clonazepam				0.5 ng/mL
- Diazepam				0.5 ng/mL
- Flunitrazepam				0.5 ng/mL
- Flurazepam				0.5 ng/mL
- Lorazepam				0.5 ng/mL
- Midazolam				0.5 ng/mL
- Nordiazepam				0.5 ng/mL
- Oxazepam				0.5 ng/mL
- Temazepam				0.5 ng/mL
- Triazolam				0.5 ng/mL
<b>Buprenorphine</b>	EIA	5 ng/mL		
- Buprenorphine				1.0 ng/mL
- Norbuprenorphine				5.0 ng/mL
<b>Cocaine</b>	EIA	20 ng/mL		
- Benzoylecgonine				4 ng/mL
- Cocaine				4 ng/mL

<b>Fentanyl</b>	EIA	5 ng/mL		
- <i>Acetylfentanyl</i>				1 ng/mL
- <i>Acrylfentanyl</i>				1 ng/mL
- <i>Butyrylfentanyl</i>				1 ng/mL
- <i>Fentanyl</i>				1 ng/mL
- <i>Furanylfentanyl</i>				1 ng/mL
<b>Marijuana (THC)</b>	EIA	4 ng/mL		1 ng/mL
<b>Methadone</b>	EIA	50 ng/mL		10 ng/mL
<b>Opiates</b>	EIA	40 ng/mL		
- <i>Codeine</i>				8 ng/mL
- <i>Morphine</i>				8 ng/mL
- <i>Hydrocodone</i>				8 ng/mL
- <i>Hydromorphone</i>				8 ng/mL
- <i>6-monoacetylmorphine (6-MAM)</i>				4 ng/mL
<b>Oxycodone</b>	EIA	40 ng/mL		
- <i>Oxycodone</i>				8 ng/mL
- <i>Oxymorphone</i>				8 ng/mL
<b>Phencyclidine (PCP)</b>	EIA	10 ng/mL		5 ng/mL
<b>Tramadol</b>	EIA	50 ng/mL		10 ng/mL